

Mass spectrometry-based proteomics

Laurent Gatto


2025-03-10

This material available under a **creative common CC-BY** license. You are free to **share** (copy and redistribute the material in any medium or format) and **adapt** (remix, transform, and build upon the material) for any purpose, even commercially.

Biological information flow

There are three main biological entities that respectively store information, act as data intermediates, and the functional units, and the figure below show how information flows between these three levels.

DNA, that lives in the nucleus of cells, is the central information storage mechanism, and encodes the blueprint of the functional units as genes. DNA is **transcribed** into **messenger RNA (mRNA)**, that re-localises outside the nucleus and is further processed into its mature *exon-only* form after removal of the non-coding *introns* sequences. Finally, the mRNA is translated by the ribosomal machinery into **proteins** directly into the endoplasmic reticulum (ER) where they are then redirected to their final destination.



In addition to the standard information workflow where DNA is transcribed into RNA that itself is translated into proteins, there is also reverse transcription, that generates (complementary) DNA




Figure 2: Sub-cellular structure of an animal cell (Source *Cell biology* on Wikipedia).

Figure 1: Information flow in biological systems (Source *Central dogma of biology* on Wikipedia).

from and RNA molecule, as well as replication of DNA (during cell division) and RNA molecules.

Why studying proteins

Proteins as the functional units in all living organisms, and they are highly dynamic. The caterpillar and the resulting butterfly have the same genome. The complement of all the expressed proteins, termed the proteome is however very different.



Figure 3: The metamorphosis from a caterpillar to a monarch butterfly. (Image from Phys.org)

There are different modalities of the proteome that are of interest. In addition to the presence or amount of protein in a biological samples, it is also important to study the interactions between proteins forming protein-protein complexes, the presence of post-transcriptional modification (such as, for example, phosphorylations), the rate at which proteins are produced and degraded, or where the proteins reside inside a cell.

The technique of choice to study proteins in a high throughput way is *mass spectrometry*.

A beginner's guide to mass spectrometry-based proteomics (Sinha and Mann 2020) is an approachable introduction to sample preparation, mass spectrometry and data analysis.

Setup

We are going to use the Bioconductor (Huber et al. 2015) MSnbase package (Laurent Gatto and Lilley 2012; Laurent Gatto, Gibb, and Rainer 2020), which can be install with the BiocManager package, available from CRAN. If BiocManager isn't available on your computer, install it with:

```
install.packages("BiocManager")
```

Now, install MSnbase and its dependencies with



Figure 4: The 'MSnbase' package.

```
BiocManager::install("MSnbase")
```

For additional information on how to analyse mass spectrometry-based proteomics data, refer to (L. Gatto and Christoforou 2014) and (Laurent Gatto 2019), or explore the the proteomics- and mass spectrometry-related packages on the Bioconductor page

The recent R for Mass Spectrometry initiative, and the documentation book provide further details and state-of-the-art infrastructure for MS and proteomics data analysis.

How does mass spectrometry work?

Mass spectrometry (MS) is a technology that *separates* charged molecules (ions) based on their mass to charge ratio (M/Z). It is often coupled to chromatography (liquid LC, but can also be gas-based GC). The time an analyte takes to elute from the chromatography column is the *retention time*.



Figure 5: The *R for Mass Spectrometry* initiative.




Figure 6: A chromatogram, illustrating the total amount of analytes over the retention time.

An mass spectrometer is composed of three components:

1. The *source*, that ionises the molecules: examples are Matrix-assisted laser desorption/ionisation (MALDI) or electrospray ionisation. (ESI)
2. The *analyser*, that separates the ions: Time of flight (TOF) or Orbitrap.
3. The *detector* that quantifies the ions.

When using mass spectrometry for proteomics, the proteins are first digested with a protease such as trypsin. In mass shotgun proteomics, the analytes assayed in the mass spectrometer are peptides.

Often, ions are subjected to more than a single MS round. After a first round of separation, the peaks in the spectra, called MS₁ spectra, represent peptides. At this stage, the only information we possess about these peptides are their retention time and their mass-to-charge (we can also infer their charge by inspecting their isotopic envelope, i.e. the peaks of the individual isotopes, see below), which is not enough to infer their identity (i.e. their sequence).

In MSMS (or MS₂), the settings of the mass spectrometer are set automatically to select a certain number of MS₁ peaks (for example 20). Once a narrow M/Z range has been selected (corresponding to one high-intensity peak, a peptide, and some background noise), it is fragmented (using for example collision-induced dissociation (CID), higher energy collisional dissociation (HCD) or electron-transfer dissociation (ETD)). The fragment ions are then themselves separated in the analyser to produce a MS₂ spectrum. The unique fragment ion pattern can then be used to infer the peptide sequence using *de novo* sequencing (when the spectrum is of high enough quality) or using a search engine such as, for example Mascot, MSGF+, . . . , that will match the observed, experimental spectrum to theoretical spectra (see details below).




Figure 7: Schematics of a mass spectrometer and two rounds of MS.

The animation below shows how 25 different ions (i.e. having different M/Z values) are separated throughout the MS analysis and are eventually detected (i.e. quantified). The final frame shows the hypothetical spectrum.

The figures below illustrate the two rounds of MS. The spectrum on the left is an MS₁ spectrum acquired after 21 minutes and 3 seconds of elution. 10 peaks, highlighted by dotted vertical lines, were selected for MS₂ analysis. The peak at M/Z 460.79 (488.8) is highlighted by a red (orange) vertical line on the MS₁ spectrum and the fragment spectra are shown on the MS₂ spectrum on the top (bot-




Figure 8: Separation and detection of ions in a mass spectrometer.

tom) right figure.

The figures below represent the 3 dimensions of MS data: a set of spectra (M/Z and intensity) of retention time, as well as the interleaved nature of MS1 and MS2 (and there could be more levels) data.




Figure 9: Parent ions in the MS1 spectrum (left) and two sected fragment ions MS2 spectra (right).




Figure 10: MS1 spectra (blue) over retention time (left). MS2 spectra (pink) interleaved between two MS1 spectra (right),

Reading and accessing MS data

Let's read a very small raw MS data file into R using the `readMSData` from the `MSnbase` package. The file that we are going to load is also available in the package.

1. Load the `MSnbase` package

```
library("MSnbase")
```

```
## Warning: multiple methods tables found for 'intersect'
```

```
## Warning: multiple methods tables found for 'intersect'
```

2. Get the path to the `dummyiTRAQ.mzXML` file

```
rawf <- dir(system.file(package = "MSnbase", dir = "extdata"),
            full.name = TRUE,
            pattern = "mzXML$")
```

```
basename(rawf)
```

```
## [1] "dummyiTRAQ.mzXML"
```

3. Read it in using the `readMSData` function.

```
x <- readMSData(rawf)
```

```
x
```

```
## MSn experiment data ("MSnExp")
```

```
## Object size in memory: 0.18 Mb
```

```
## - - - Spectra data - - -
```

```
## MS level(s): 2
```

```
## Number of spectra: 5
```

```
## MSn retention times: 25:01 - 25:02 minutes
```

```
## - - - Processing information - - -
```

```
## Data loaded: Mon Mar 10 11:47:33 2025
```

```
## MSnbase version: 2.32.0
```

```
## - - - Meta data - - -
```

```
## phenoData
```

```
##   rowNames: dummyiTRAQ.mzXML
```

```
##   varLabels: sampleNames
```

```
##   varMetadata: labelDescription
```

```
## Loaded from:
```

```
##   dummyiTRAQ.mzXML
```

```
## protocolData: none
```

```
## featureData
```

```
##   featureNames: F1.S1 F1.S2 ... F1.S5 (5 total)
```

```
##   fvarLabels: spectrum
```

```
##   fvarMetadata: labelDescription
```

```
## experimentData: use 'experimentData(object)'
```

The object that is returned by `readMSData` is of class `MSnExp`, that can store, access and manipulate raw MS data. Note that here we are focusing on MS-based proteomics data, but this also applied to MS-based metabolomics data.

```
class(x)
```

```
## [1] "MSnExp"
## attr(,"package")
## [1] "MSnbase"
```

4. We can find out how many spectra are available in that data using the function `length`. Full MS acquisitions would contain hundreds of thousands spectra.

```
length(x)
```

```
## [1] 5
```

5. We can use various accessor function to get the MS level of these spectra, their retention time, or the M/Z and intensity of the precursor peaks of the ion corresponding to the MS2 spectra.

```
msLevel(x)
```

```
## F1.S1 F1.S2 F1.S3 F1.S4 F1.S5
##      2      2      2      2      2
```

```
rtime(x)
```

```
## F1.S1 F1.S2 F1.S3 F1.S4 F1.S5
## 1501.35 1501.59 1501.85 1502.07 1502.31
```

```
precursorMz(x)
```

```
## F1.S1 F1.S2 F1.S3 F1.S4 F1.S5
## 645.3741 546.9586 645.3741 716.3405 437.8040
```

```
precursorIntensity(x)
```

```
## F1.S1 F1.S2 F1.S3 F1.S4 F1.S5
## 47659400 26356100 23432400 24854800 7052960
```

6. We can also extract individual spectra using `[]` and plot them.

```
x[[3]]
```

```
## Object of class "Spectrum2"
## Precursor: 645.3741
## Retention time: 25:02
## Charge: 2
## MSn level: 2
## Peaks count: 2125
## Total ion count: 150838188
```



```
plot(x[[3]])
```




Figure 11: Visualisation of the 3rd MS spectrum in our small test data set.

Exercise

For the rest of this tutorial, we will be using a slightly larger dataset (still tiny compared to full acquisitions) that is distributed with the MSnbase package. Load it as shown below and compute the number of spectra available in that dataset, their MS level, and the retention time range over which these spectra have been acquired.

```
data(itraqdata)

length(itraqdata)

## [1] 55

unique(msLevel(itraqdata))

## [1] 2

formatRt(range(rtime(itraqdata)))

## [1] "19:09" "50:18"
```

This object also contains additional metadata for each spectrum, that can be accessed, as a `data.frame`, with `fData`.

Identification

The raw data is still a long way of obtaining biologically relevant proteomics data. The first step to obtain proteomics data is to identify the peptides that have been acquired in the MS. Peptide identification work by comparing expected and observed spectra. As shown below, when a precursor peptide ion is fragmented in a CID cell, it breaks at specific bonds, producing sets of peaks (a , b , c and x , y , z) that can be predicted.

It is thus possible to calculate the expected set of fragment peaks for a given peptide, such as *SIGFEGDSIGR* below.

```
calculateFragments("SIGFEGDSIGR")
```

```
##           mz ion type pos z      seq
## 1  88.03931 b1    b   1 1        S
## 2 201.12337 b2    b   2 1       SI
## 3 258.14483 b3    b   3 1      SIG
## 4 405.21324 b4    b   4 1     SIGF
## 5 534.25583 b5    b   5 1    SIGFE
## 6 591.27729 b6    b   6 1   SIGFEG
## 7 706.30423 b7    b   7 1  SIGFEGD
## 8 793.33626 b8    b   8 1 SIGFEGDS
## 9 906.42032 b9    b   9 1 SIGFEGDSI
## 10 963.44178 b10   b  10 1 SIGFEGDSIG
## 11 175.11895 y1    y   1 1        R
## 12 232.14041 y2    y   2 1       GR
## 13 345.22447 y3    y   3 1      IGR
## 14 432.25650 y4    y   4 1     SIGR
## 15 547.28344 y5    y   5 1    DSIGR
## 16 604.30490 y6    y   6 1   GDSIGR
## [ reached 'max' / getOption("max.print") -- omitted 16 rows ]
```

The last step is to compare observed and expected peaks. If there is a good match, the MS2 spectrum is assigned the peptide sequence.


```
itraqdata2 <- pickPeaks(itraqdata, verbose = FALSE)
s <- "SIGFEGDSIGR"
plot(itraqdata2[[14]], s, main = s)
```

It is also possible to plot 2 spectra to compare them directly.

```
plot(itraqdata2[[25]], itraqdata2[[28]],
     sequences = rep("IMIDLDGTENK", 2))
```

In a full experiment, all possible peptides from the known (or relevant) proteome of interest (such as databases that can be downloaded

Figure 12: Peptide fragmentation.



Biemann, K *Methods Enzymol* (1990) **193** 886-8




Figure 13: Matching observed and expected peaks.




Figure 14: Direct comparison of 2 MS2 spectra.

from the UniProt site¹) are compared to the millions of observed spectra.


¹ The Universal Protein Resource (UniProt) is a freely and accessible comprehensive resource for protein sequence and annotation data.

From the list of identified peptides, it is then necessary to infer the most provable proteins that were present in the biological sample.

Exercise

Plot the 44th spectrum of the itraqdata2 experiment. The sequence can be accessed in the feature metadata with

```
fData(itraqdata2)$PeptideSequence[[44]]  
  
## [1] LDVLDGLK  
## 47 Levels: AADALLLK AAGHDGK AAKPEAPAASPAPALGAR AALESAPK ... VVVGSK  
  
plot(itraqdata2[[44]], s,  
     main = fData(itraqdata2)$PeptideSequence[[44]])
```



Quantitation


The last step of MS data processing is to quantify peptide abundances in the biological samples. The table below summarises the different possibilities depending whether the proteins or peptides are labelled, and whether the quantitation is performed in MS₁ or MS₂.

	Label-free	Labelled
MS ₁	XIC	SILAC, ¹⁵ N
MS ₂	Counting	iTRAQ, TMT

Label-free MS₂: Spectral counting


In spectral counting, one simply counts the number of quantified peptides that are assigned to a protein.

Statistical challenge The matching between millions of observed and possible spectra causes real challenges due to the large search space and the risk of false positives. See [Kall:2008] for further reading.
Statistical challenge Protein inference is a difficult task, as peptides often match multiple proteins (either different isoforms or proteins stemming from different genes but with identical domains), which leads to the definition of protein groups, i.e. sets of proteins that can't be distinguished with the set of identified peptides at hand. See [Kozlowski:2005] for further reading.



Labelled MS2: Isobaric tagging

Isobaric tagging refers to the labelling using isobaric tags, i.e. chemical tags that have the same mass and hence can't be distinguished by the spectrometer. The peptides of different samples (4, 6, 10 or 11) are labelled with different tags and combined prior to mass spectrometry acquisition. Given that they are isobaric, all identical peptides, irrespective of the tag and thus the sample of origin, are co-analysed, up to fragmentation prior to MS2 analysis. During fragmentation, the isobaric tags fall off, fragment themselves, and result in a set of sample-specific peaks. These specific peaks can be used to infer sample-specific quantitation, while the rest of the MS2 spectrum is used for identification.



Label-free MS1: extracted ion chromatograms

In label-free quantitation, the precursor peaks that match an identified peptide are integrated of retention time and the area under that *extracted ion chromatogram* is used to quantify that peptide in that sample.




Figure: credit Johannes Rainer.

Labelled MS1: SILAC

In SILAC quantitation, sample are grown in a medium that contains heavy amino acids (typically arginine and lysine). All proteins grown

in this *heavy* growth medium contain the heavy form of these amino acids. Two samples, one grown in heavy medium, and one grown in normal (light) medium are then combined and analysed together. The heavy peptides precursor peaks are systematically shifted compared to the light ones, and the ratio between the height of a heavy and light peaks can be used to calculate peptide and protein fold-changes.




Figure: credit Wikimedia Commons.

Exercise

As its name implies, the `itraqdata` is an iTRAQ-based isobar quantitation experiment. We can visualise the reporter peaks as follows:

```
plot(itraqdata[[14]], reporters = iTRAQ4, full = TRUE)
```




Figure 15: Visualisation of the iTRAQ reporter peaks.

Interpret the figure above.

From the closed-up panel, we see that the 4 reporter ions indicate that there was slightly less of that peptide in the two first samples, i.e. those labelled with the 114 and 115 reporters.

The rest of the spectrum (i.e $m/z > 120$) correspond to the dissociated fragment ions, used for identification.

Quantification

We can quantify these four peaks with the `quantify` method, to produce an object of class `MSnSet` containing quantitation data. The quantitation values can be accessed with `exprs`. This data also contains feature metadata that can be accessed with the `fData` function.

```
msnset <- quantify(itraqdata, method = "trap",
                  reporters = iTRAQ4)

msnset

## MSnSet (storageMode: lockedEnvironment)
## assayData: 55 features, 4 samples
##   element names: exprs
## protocolData: none
## phenoData
##   sampleNames: iTRAQ4.114 iTRAQ4.115 iTRAQ4.116 iTRAQ4.117
##   varLabels: mz reporters
##   varMetadata: labelDescription
## featureData
##   featureNames: X1 X10 ... X9 (55 total)
##   fvarLabels: spectrum ProteinAccession ... collision.energy (15 total)
##   fvarMetadata: labelDescription
## experimentData: use 'experimentData(object)'
## Annotation: No annotation
## - - - Processing information - - -
## Data loaded: Wed May 11 18:54:39 2011
## Updated from version 0.3.0 to 0.3.1 [Fri Jul  8 20:23:25 2016]
## iTRAQ4 quantification by trapezoidation: Mon Mar 10 11:47:39 2025
## MSnbase version: 1.1.22

head(exprs(msnset))

##      iTRAQ4.114 iTRAQ4.115 iTRAQ4.116 iTRAQ4.117
## X1    1347.6158  2247.3097  3927.6931  7661.1463
## X10    739.9861   799.3501   712.5983   940.6793
## X11  27638.3582 33394.0252 32104.2879 26628.7278
## X12  31892.8928 33634.6980 37674.7272 37227.7119
## X13  26143.7542 29677.4781 29089.0593 27902.5608
## X14   6448.0829  6234.1957  6902.8903  6437.2303
```

```
head(fData(msnset))
```

```
##      spectrum ProteinAccession      ProteinDescription
## X1      1      BSA      bovine serum albumin
## X10     10     ECA1422 glucose-1-phosphate cytidylyltransferase
## X11     11     ECA4030      50S ribosomal subunit protein L4
## X12     12     ECA3882      chaperone protein DnaK
## X13     13     ECA1364      succinyl-CoA synthetase alpha chain
## X14     14     ECA0871      NADP-dependent malic enzyme
##      PeptideSequence fileIdx retention.time precursor.mz precursor.intensity
## X1      NYQEAK      1      1149.31      520.7833      3449020
## X10     VTLVDTGEHSMGTGGR      1      1503.03      573.9539      7849420
## X11      SPIWR      1      1663.61      401.7392      41253600
## X12     TAIDDALK      1      1663.86      567.8339      23549500
## X13      SILINK      1      1664.08      488.3269      13025200
## X14     DFEVVNNESDPR      1      1664.32      782.8715      18405000
##      charge peaks.count      tic ionCount ms.level acquisition.number
## X1      2      1922 26413754 26413754      2      2
## X10     3      1376 24482281 24482281      2      11
## X11     2      1571 231075934 231075934      2      12
## X12     2      2397 247323187 247323187      2      13
## X13     2      2574 207247502 207247502      2      14
## X14     2      1829 115317275 115317275      2      15
##      collision.energy
## X1      40
## X10     40
## X11     40
## X12     40
## X13     40
## X14     40
```

Quantitative data processing

In our examples, we not have processing data for the 55 peptides and 4 samples. In this data, there is only 1 missing value, corresponding to an absent reporter peak. We are going to simply drop that feature.

```
table(is.na(exprs(msnset)))
```

```
##
## FALSE TRUE
## 219    1
```

```
msnset <- filterNA(msnset)
```

Statistical challenge Quantitative data processing come with numerous statistical challenges such as missing data handling, aggregation of quantitation data, data normalisation, ... and the effect of these on the downstream statistical tests.

In MS1 label-free experiments, given that each sample is acquired independently, the proportion of missing values can be as high several tens of percent. In such situations, removing rows with missing values isn't possible at all. Imputation is possible, albeit tricky, as different mechanisms can be responsible for missing value that appear either at random or not at random (Lazar et al. 2016).

Next, we aggregate the spectrum-level quantitation values into protein-level data using the median and the `combineFeatures` function:

```
prots <- combineFeatures(msnset, fcol = "ProteinAccession",
                        method = "median")
head(exprs(prots))
```

##	iTRAQ4.114	iTRAQ4.115	iTRAQ4.116	iTRAQ4.117
## BSA	1347.616	2247.310	3927.693	7661.1463
## ECA0172	17593.548	18545.620	19361.837	18328.2365
## ECA0435	4923.628	5557.818	5775.203	5079.2952
## ECA0452	1524.148	1399.897	1547.218	1563.2299
## ECA0469	1069.945	1035.689	1029.420	999.6957
## ECA0621	1101.062	1124.167	1140.093	1191.8055

Following on from here, many data processing such as normalisation, non-specific filtering, and hypothesis testing is very similar to other omics data.

Applications in statistical learning

In this section, I illustrate a use case of mass spectrometry-based proteomics to infer sub-cellular protein localisation and the application of statistical learning. This section requires the `pRoloc` and `pRolocdata` packages (L. Gatto et al. 2014; Breckels et al. 2016). Both packages can be installed with the `BiocManager::install` function, as illustrated above.

```
library("pRoloc")

## Warning: multiple methods tables found for 'intersect'

library("pRolocdata")
```

The hyperLOPIT2015 data contains localisation data for 5032 proteins in mouse embryonic stem cells (Christoforou et al. 2016). The protein information is collected along a 20 fraction density gradient - our data matrix has dimensions 5032 rows and 20 columns. We use PCA to easily visualise it in 2 dimensions.




Figure 16: The 'pRoloc' package.

```
data(hyperLOPIT2015)
## set colours
setStockcol(paste0(getStockcol(), 80))
## produce PCA plot
plot2D(hyperLOPIT2015)
```




Figure 17: Mouse stem cell spatial proteomics data from Christoforou et al.

Each dot on the PCA plot represents a single protein. The protein is either of unknown localisation, and represented by a grey circle, or is of known localisation (and called an organelle marker) and is coloured according to its expected sub-cellular localisation.

In the next figure, we have trained a support vector machine (SVM) model and classified proteins of unknown localisation using an SVM model. The size of the dots are scaled according to the classifier score and only assignments corresponding to a false discovery rate of 5% have been considered.

```
sz <- exp(fData(hyperLOPIT2015)$svm.score) - 1
plot2D(hyperLOPIT2015, fcol = "final.assignment", cex = sz)
```

Session information

```
## R version 4.4.1 (2024-06-14)
```




Figure 18: New sub-cellular assignment after using support vector machine classifier.

```
## Platform: x86_64-pc-linux-gnu
## Running under: Ubuntu 24.04.2 LTS
##
## Matrix products: default
## BLAS: /opt/Rpatched/lib/R/lib/libRblas.so
## LAPACK: /opt/Rpatched/lib/R/lib/libRlapack.so; LAPACK version 3.12.0
##
## locale:
## [1] LC_CTYPE=en_US.UTF-8 LC_NUMERIC=C
## [3] LC_TIME=en_US.UTF-8 LC_COLLATE=en_US.UTF-8
## [5] LC_MONETARY=en_US.UTF-8 LC_MESSAGES=en_US.UTF-8
## [7] LC_PAPER=en_US.UTF-8 LC_NAME=C
## [9] LC_ADDRESS=C LC_TELEPHONE=C
## [11] LC_MEASUREMENT=en_US.UTF-8 LC_IDENTIFICATION=C
##
## time zone: Europe/Brussels
## tzcode source: system (glibc)
##
## attached base packages:
## [1] stats4 stats graphics grDevices utils datasets methods
## [8] base
##
```

```

## other attached packages:
## [1] pRolocdata_1.44.1      pRoloc_1.46.1      BiocParallel_1.40.0
## [4] MLInterfaces_1.86.0    cluster_2.1.6      annotate_1.84.0
## [7] XML_3.99-0.18          AnnotationDbi_1.68.0 IRanges_2.40.1
## [10] MSnbase_2.32.0         ProtGenerics_1.38.0 S4Vectors_0.44.0
## [13] mzR_2.40.0             Rcpp_1.0.14        Biobase_2.66.0
## [16] BiocGenerics_0.52.0
##
## loaded via a namespace (and not attached):
## [1] splines_4.4.1          filelock_1.0.3
## [3] tibble_3.2.1           hardhat_1.4.1
## [5] preprocessCore_1.68.0  pROC_1.18.5
## [7] rpart_4.1.23           lifecycle_1.0.4
## [9] httr2_1.1.0            doParallel_1.0.17
## [11] globals_0.16.3        lattice_0.22-6
## [13] MASS_7.3-60.2          MultiAssayExperiment_1.32.0
## [15] dendextend_1.19.0      magrittr_2.0.3
## [17] limma_3.62.2           plotly_4.10.4
## [19] rmarkdown_2.29         yaml_2.3.10
## [21] MsCoreUtils_1.18.0     DBI_1.2.3
## [23] RColorBrewer_1.1-3     lubridate_1.9.4
## [25] abind_1.4-8            zlibbioc_1.52.0
## [27] GenomicRanges_1.58.0   purrr_1.0.4
## [29] mixtools_2.0.0         AnnotationFilter_1.30.0
## [31] nnet_7.3-20            rappdirs_0.3.3
## [33] ipred_0.9-15           lava_1.8.1
## [35] GenomeInfoDbData_1.2.13 listenv_0.9.1
## [37] parallelly_1.42.0      ncdf4_1.23
## [39] codetools_0.2-20       DelayedArray_0.32.0
## [41] xml2_1.3.7             tidyselect_1.2.1
## [43] UCSC.utils_1.2.0       farver_2.1.2
## [45] viridis_0.6.5          matrixStats_1.5.0
## [47] BiocFileCache_2.14.0   jsonlite_1.9.1
## [49] caret_7.0-1            e1071_1.7-16
## [51] survival_3.6-4         iterators_1.0.14
## [53] foreach_1.5.2          segmented_2.1-4
## [55] tools_4.4.1            progress_1.2.3
## [57] glue_1.8.0             prodlim_2024.06.25
## [59] gridExtra_2.3          SparseArray_1.6.2
## [61] xfun_0.51              tuftes_0.13
## [63] MatrixGenerics_1.18.1  GenomeInfoDb_1.42.3
## [65] dplyr_1.1.4            withr_3.0.2
## [67] BiocManager_1.30.25    fastmap_1.2.0
## [69] digest_0.6.37          timechange_0.3.0

```

```
## [71] R6_2.6.1                colorspace_2.1-1
## [73] gtools_3.9.5             lpSolve_5.6.23
## [75] biomaRt_2.62.1           RSQLite_2.3.9
## [77] tidyr_1.3.1              generics_0.1.3
## [79] hexbin_1.28.5            data.table_1.17.0
## [81] recipes_1.1.1            FNN_1.1.4.1
## [83] class_7.3-22             prettyunits_1.2.0
## [85] PSMatch_1.10.0           httr_1.4.7
## [87] htmlwidgets_1.6.4        S4Arrays_1.6.0
## [89] ModelMetrics_1.2.2.2     pkgconfig_2.0.3
## [91] gtable_0.3.6             timeDate_4041.110
## [93] blob_1.2.4               impute_1.80.0
## [95] XVector_0.46.0           htmltools_0.5.8.1
## [97] MALDIquant_1.22.3        clue_0.3-66
## [99] scales_1.3.0             png_0.1-8
## [ reached getOption("max.print") -- omitted 52 entries ]
```

References

- Breckels, L M, C M Mulvey, K S Lilley, and L Gatto. 2016. "A Bioconductor Workflow for Processing and Analysing Spatial Proteomics Data." *F1000Res* 5: 2926. <https://doi.org/10.12688/f1000research.10411.2>.
- Christoforou, A, C M Mulvey, L M Breckels, A Geladaki, T Hurrell, P C Hayward, T Naake, et al. 2016. "A Draft Map of the Mouse Pluripotent Stem Cell Spatial Proteome." *Nat Commun* 7 (January): 8992. <https://doi.org/10.1038/ncomms9992>.
- Gatto, Laurent. 2019. *Bioconductor Tools for Mass Spectrometry and Proteomics*. <https://rawgit.com/lgatto/bioc-ms-prot/master/lab.html>.
- Gatto, Laurent, Sebastian Gibb, and Johannes Rainer. 2020. "MSnbase, Efficient and Elegant R-Based Processing and Visualisation of Raw Mass Spectrometry Data." *bioRxiv*. <https://doi.org/10.1101/2020.04.29.067868>.
- Gatto, Laurent, and Kathryn S Lilley. 2012. "MSnbase-an R/Bioconductor Package for Isobaric Tagged Mass Spectrometry Data Visualization, Processing and Quantitation." *Bioinformatics* 28 (2): 288–89.
- Gatto, L, L M Breckels, S Wiczorek, T Burger, and K S Lilley. 2014. "Mass-Spectrometry-Based Spatial Proteomics Data Analysis Using pRoloc and pRolocdata." *Bioinformatics* 30 (9): 1322–24. <https://doi.org/10.1093/bioinformatics/btu013>.
- Gatto, L, and A Christoforou. 2014. "Using R and Bioconductor for Proteomics Data Analysis." *Biochim. Biophys. Acta* 1844 (1 Pt A): 42–51.

- Huber, W, V J Carey, R Gentleman, S Anders, M Carlson, B S Carvalho, H C Bravo, et al. 2015. "Orchestrating High-Throughput Genomic Analysis with Bioconductor." *Nat. Methods* 12 (2): 115–21.
- Lazar, C, L Gatto, M Ferro, C Bruley, and T Burger. 2016. "Accounting for the Multiple Natures of Missing Values in Label-Free Quantitative Proteomics Data Sets to Compare Imputation Strategies." *J Proteome Res* 15 (4): 1116–25. <https://doi.org/10.1021/acs.jproteome.5b00981>.
- Sinha, Ankit, and Matthias Mann. 2020. "A beginner's guide to mass spectrometry-based proteomics." *The Biochemist*, September. <https://doi.org/10.1042/BI020200057>.